

### **Sporulation plates**

Potassium acetate	1 % w/v
Bacto Agar	2 % w/v
Amino acid stock (sporulation plate)	0.25 % v/v

### **Amino acid stock (sporulation plate)**

–dissolved in 100 mL H<sub>2</sub>O

Uracil, Histidine, Methionine	0.2 g
Leucine, Lysine	0.3 g

### **Antibiotic selection plates**

Hygromycin B (Hph) (InvivoGen)	300 mg/L
Geneticin disulphate salt (G418) (Sigma)	200 mg/L
Nourseothricin (clonNAT) (Werner BioAgents)	100 mg/L